

## Home Bed Assist Handle

HCPCS:




FIGURE A



FIGURE B

- For use with most home mattress/box spring/bed frame type beds
- Weight supported by floor on rubber tip feet, not just by box spring and mattress
- Adjustable handle height and leg extensions to accommodate a wide range of bed heights (see specifications for exact fit)
- Simple to use handle assembly slides under mattress to provide a safe and secure grip for assistance getting into and out of bed
- Non-Slip foam cover on horizontal stabilizer bars provides increased stability within the mattress
- Folds flat for easy storage and portability

### PRODUCTS

### Description

RTL15063-ADJ	Bed Assist Rail,Adj.Height,RTL,1 EA
--------------	-------------------------------------

### SPECIFICATIONS

### Description

Height Adjustment	15" - 20"
Rail Dimensions	32.75" (D) x 19.75" (W) x 13.5-16.5" (H)
Warranty	Limited Lifetime
Weight Capacity	250 lbs