

Transfer Tub Bench

HCPCS: E0248,E0248



FIGURE A



FIGURE B



FIGURE C

- Tool-free assembly of back, legs and arm (Figure A - C)
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Height adjusts in $\frac{1}{2}$ " increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

PRODUCTS

Description

12011KD-1	Transfer Bench KD,Tool Free, 1/cs
-----------	-----------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	13 lbs. 13 lbs. 24 lbs.
HCPCS	E0248
Outside Legs	32" (W) x 23" (D)
Seat Dimensions	26" (W) x 18.5" - 19.5" (D)
Seat Height	17.5" - 22.5"
Warranty	Limited Lifetime
Weight	10 lbs.
Weight Capacity	400 lbs

ACCESSORIES

Description

12011SS	SuctionCups,Sm.12011,1pr/bx, BX
12011SL	SuctionCups,Lg.12011,1pr/bx, BX