

## MODEL

Jazzy Elite HD-2HD-SS	K0824
Jazzy Elite HD-2HD-C	K0825
<b>Weight Capacity</b>	450 lbs.
<b>Turning Radius</b>	24"
<b>Overall Width<sup>2</sup></b>	24.25"
<b>Overall Length<sup>2</sup></b>	35.5"
<b>Maximum Speed<sup>1</sup></b>	4 mph
<b>Ground Clearance<sup>2</sup></b>	3"
<b>Front Wheels</b>	3" solid anti-tips
<b>Drive Wheels</b>	14" solid
<b>Rear Wheels</b>	6" solid casters
<b>Drivetrain</b>	Two motor, front-wheel drive
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Suspension Type</b>	Limited
<b>High-back Seat Sizes</b>	18" x 18-20", 20" x 18-20" and 22" x 20-22"
<b>Seat-to-Floor Range</b>	20"-22"
<b>Per Charge Range<sup>2,4</sup></b>	Up to 9.8 miles
<b>Standard Electronics</b>	60A, Dynamic Shark controller
<b>Component Weights</b>	Base 99 lbs. (heaviest piece) Seat 52 lbs. Batteries 37 lbs. each
<b>Battery Requirements<sup>3,5</sup></b>	(2) 12 volt, NF-22 batteries
<b>Battery Charger</b>	5 amp, off-board
<b>Warranty</b>	Frame: Lifetime limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>5</sup> Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

**PRIDE**  
MOBILITY  
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

JAZZYELITEHD\_060719

**Jazzy**  
ELITE HD

**HEAVY DUTY  
FRONT-WHEEL  
DRIVE**



**AMERICA'S FAVORITE POWER CHAIR**

gojazzy.com |    

# ENHANCED OUTDOOR PERFORMANCE

High-back seat

Intuitive controller  
is easy to operate

Color-through shroud

Easy rear  
access to  
batteries

Angle-adjustable  
foot platform offers  
greater comfort

Large 14" drive wheels



Front-wheel drive design




Dynamic Shark Controller

The **Jazzy® Elite HD** features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front anti-tip wheels for handling various terrains with ease.

- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 70" lap belt

## Available Colors



Jazzy Red

Jazzy Blue

## Optional Seating

Deluxe contour, limited recline, high back with solid seat pan, available sizes:

- 18" x 18"
- 18" x 20"

## Options

- Bell shaped foot platform
- Elevating leg rests
- Heel loops
- Height adjustable foot platform bracket
- Residual limb support
- Swing-away legrests
- Swing-away joystick

## Accessories

- Cane/crutch holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger